



# PUB3 Antibody

|                     |  |
|---------------------|--|
| Product Code        | CSB-PA816753XA01DOA  |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| Uniprot No.         | Q8GWV5   |
| Immunogen           | Recombinant Arabidopsis thaliana PUB3 protein                                    |
| Raised In           | Rabbit   |
| Species Reactivity  | Arabidopsis thaliana   |
| Tested Applications | ELISA, WB (ensure identification of antigen)                                     |
| Form                | Liquid   |
| Conjugate           | Non-conjugated   |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified  |
| Isotype             | IgG  |
| Clonality           | Polyclonal   |
| Product Type        | Polyclonal Antibody  |
| Immunogen Species   | Arabidopsis thaliana (Mouse-ear cress)   |
| Lead Time           | Made-to-order (14-16 weeks)  |
| Gene Names          | PUB3   |