



Recombinant *Geobacillus stearothermophilus* Glyceraldehyde-3-phosphate dehydrogenase (gap)

| | |
|--------------------------|---|
| Product Code | CSB-EP009232GFM-B |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P00362 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Geobacillus stearothermophilus</i> (<i>Bacillus stearothermophilus</i>) |
| Purity | >85% (SDS-PAGE) |
| Sequence | AVKVGINGF GRIGRNVFRA ALKNPDIEVV AVNDLTDANT LAHLLKYDSV HGRLDAEVS SV NGNNLVVNGK EIIVKAERDP ENLAWGEIGV DIVVESTGRF TKREDAAKHL EAGAKKVIIS APAKNEDITI VMGVNQDKYD PKAHHVISNA SCTTNCLAPF AKVLHEQFGI VRGMMTTVHS YTNDQRILD L PHKDLRRARA AAESIIP TTT GAAKAVALVL PELKGKLN GM AMRVPTPNVS VVDLVAELEK EVTVEEVNAA LKAAAEGELK GILAYSEEPL VSRDYNGSTV SSTIDALSTM VIDGKMVKVV SWYDNETGYS HRVVDLAAYI ASKGL |
| Source | E.coli |
| Target Names | gap |
| Protein Names | Recommended name: Glyceraldehyde-3-phosphate dehydrogenase Short name= GAPDH EC= 1.2.1.12 |
| Expression Region | 2-335 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |