



RPL8C Antibody

Product Code	CSB-PA677072XA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q42064
Immunogen	Recombinant Arabidopsis thaliana RPL8C protein
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Lead Time	Made-to-order (14-16 weeks)
Target Names	RPL8C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.