




SBT6.1 Antibody

Product Code	CSB-PA877556XA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0WUG6
Immunogen	Recombinant Arabidopsis thaliana (Mouse-ear cress) SBT6.1 protein (175-473aa)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana (Mouse-ear cress)
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Target Names	SBT6.1 Antibody

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes: SBT6.1 Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 36.6 kDa

Observed band size: 37 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.