



SPL13A Antibody

Product Code	CSB-PA495133XA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B9DI20
Immunogen	Recombinant Arabidopsis thaliana SPL13A protein
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Lead Time	Made-to-order (14-16 weeks)
Target Names	SPL13A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.