



Y48B6A.1 Antibody

Product Code	CSB-PA880028XA01CXY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9U2A9
Immunogen	Recombinant Caenorhabditis elegans Y48B6A.1 protein
Raised In	Rabbit
Species Reactivity	Caenorhabditis elegans
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Caenorhabditis elegans
Lead Time	Made-to-order (14-16 weeks)
Target Names	Y48B6A.1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.