



Y92H12BL.1 Antibody

Product Code	CSB-PA839767XA01CXY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8MXQ7
Immunogen	Recombinant Caenorhabditis elegans Y92H12BL.1 protein
Raised In	Rabbit
Species Reactivity	Caenorhabditis elegans
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Caenorhabditis elegans
Lead Time	Made-to-order (14-16 weeks)
Target Names	Y92H12BL.1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.