



emb-27 Antibody

| | |
|----------------------------|----------------------------------------------------------------------------------|
| Product Code | CSB-PA917108XA01CXY |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | G5EG38 |
| Immunogen | Recombinant Caenorhabditis elegans emb-27 protein |
| Raised In | Rabbit |
| Species Reactivity | Caenorhabditis elegans |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Caenorhabditis elegans |
| Lead Time | Made-to-order (14-16 weeks) |
| Target Names | emb-27 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |