



CDK5RAP1 Antibody

Product Code	CSB-PA822305XA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96SZ6
Immunogen	Recombinant Homo sapiens (Human) CDK5RAP1 protein
Raised In	Rabbit
Species Reactivity	Homo sapiens (Human)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (14-16 weeks)
Target Names	CDK5RAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.