



3A3 Antibody

Product Code	CSB-PA321432XA01HSW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16253
Immunogen	Recombinant Haemonchus contortus (Barber pole worm) 3A3 protein
Raised In	Rabbit
Species Reactivity	Haemonchus contortus (Barber pole worm)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Haemonchus contortus (Barber pole worm)
Lead Time	Made-to-order (14-16 weeks)
Target Names	3A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.