



Gdpgp1 Antibody

Product Code	CSB-PA659698XA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3TLS3
Immunogen	Recombinant Mus musculus (Mouse) Gdpgp1 protein
Raised In	Rabbit
Species Reactivity	Mus musculus (Mouse)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Gdpgp1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.