



Mab21I2 Antibody

Product Code	CSB-PA804385XA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8BPP1
Immunogen	Recombinant Mus musculus (Mouse) Mab21I2 protein
Raised In	Rabbit
Species Reactivity	Mus musculus (Mouse)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Mab21I2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.