



CBG14305 Antibody

Product Code	CSB-PA228839XA01CXX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A8XJQ6
Immunogen	Recombinant Caenorhabditis briggsae CBG14305 protein
Raised In	Rabbit
Species Reactivity	Caenorhabditis briggsae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Caenorhabditis briggsae
Lead Time	Made-to-order (14-16 weeks)
Target Names	CBG14305
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.