




nagD Antibody

Product Code	CSB-PA640235XA01FLF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2FZX0
Immunogen	Recombinant Staphylococcus aureus (strain NCTC 8325) nagD protein (1-259aa)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus (strain NCTC 8325)
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Staphylococcus aureus (strain NCTC 8325)
Target Names	nagD

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes:nagD Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/20000 dilution

Predicted band size: 34.9 kDa

Observed band size: 35 kDa

Exposure time: 10s.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.