



nd2 Antibody

Product Code	CSB-PA811519XA01CXX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8HEC1
Immunogen	Recombinant Caenorhabditis briggsae nd2 protein
Raised In	Rabbit
Species Reactivity	Caenorhabditis briggsae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Caenorhabditis briggsae
Lead Time	Made-to-order (14-16 weeks)
Target Names	nd2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.