



ND4L Antibody

Product Code	CSB-PA015080XA01DOS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24874
Immunogen	Recombinant Ascaris suum (Pig roundworm) (Ascaris lumbricoides) ND4L protein
Raised In	Rabbit
Species Reactivity	Ascaris suum (Pig roundworm) (Ascaris lumbricoides)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Ascaris suum (Pig roundworm) (Ascaris lumbricoides)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ND4L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.