




pak-1 Antibody

Product Code	CSB-PA619143XA01CXY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q17850
Immunogen	Recombinant Caenorhabditis elegans pak-1 protein (293-572aa)
Raised In	Rabbit
Species Reactivity	Caenorhabditis elegans
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Caenorhabditis elegans
Target Names	pak-1

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes: pak-1 Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/20000 dilution

Predicted band size: 38.3 kDa

Observed band size: 38 kDa

Exposure time: 10s.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.