



Y54F10AM.8 Antibody

Product Code	CSB-PA860627XA01CXY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BL07
Immunogen	Recombinant Caenorhabditis elegans Y54F10AM.8 protein
Raised In	Rabbit
Species Reactivity	Caenorhabditis elegans
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Caenorhabditis elegans
Lead Time	Made-to-order (14-16 weeks)
Target Names	Y54F10AM.8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.