



Scnn1a Antibody

Product Code	CSB-PA020848XA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37089
Immunogen	Recombinant Rattus norvegicus (Rat) Scnn1a protein
Raised In	Rabbit
Species Reactivity	Rattus norvegicus (Rat)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Scnn1a
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.