



TFIP11 Antibody

Product Code	CSB-PA023435XA01FQF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A1XD97
Immunogen	Recombinant Canis lupus familiaris (Dog) (Canis familiaris) TFIP11 protein
Raised In	Rabbit
Species Reactivity	Canis lupus familiaris (Dog) (Canis familiaris)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Lead Time	Made-to-order (14-16 weeks)
Target Names	TFIP11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.