



tfip11 Antibody

Product Code	CSB-PA606592XA01XBF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0IIX9
Immunogen	Recombinant Xenopus tropicalis (Western clawed frog) (Silurana tropicalis) tfip11 protein
Raised In	Rabbit
Species Reactivity	Xenopus tropicalis (Western clawed frog) (Silurana tropicalis)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Xenopus tropicalis (Western clawed frog) (Silurana tropicalis)
Lead Time	Made-to-order (14-16 weeks)
Target Names	tfip11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.