



tcaA Antibody

Product Code	CSB-PA748560XA01SKY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7A3X6
Immunogen	Recombinant Staphylococcus aureus (strain N315) tcaA protein
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus (strain N315)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Staphylococcus aureus (strain N315)
Lead Time	Made-to-order (14-16 weeks)
Target Names	tcaA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.