



TFIP11 Antibody

Product Code	CSB-PA023435XA01EQV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A1XD93
Immunogen	Recombinant Pan troglodytes (Chimpanzee) TFIP11 protein
Raised In	Rabbit
Species Reactivity	Pan troglodytes (Chimpanzee)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pan troglodytes (Chimpanzee)
Lead Time	Made-to-order (14-16 weeks)
Target Names	TFIP11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.