



Tfip11 Antibody

Product Code	CSB-PA713160XA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5U2Y6
Immunogen	Recombinant Rattus norvegicus (Rat) Tfip11 protein
Raised In	Rabbit
Species Reactivity	Rattus norvegicus (Rat)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Tfip11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.