



Vesicular acetylcholine transporter Antibody

Product Code	CSB-PA841966XA01FQG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q91514
Immunogen	Recombinant Torpedo torpedo (Common torpedo) (Torpedo ocellata) Vesicular acetylcholine transporter protein
Raised In	Rabbit
Species Reactivity	Torpedo torpedo (Common torpedo) (Torpedo ocellata)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Torpedo torpedo (Common torpedo) (Torpedo ocellata)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.