



ZK1290.11 Antibody

Product Code	CSB-PA788806XA01CXY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q09338
Immunogen	Recombinant Caenorhabditis elegans ZK1290.11 protein
Raised In	Rabbit
Species Reactivity	Caenorhabditis elegans
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Caenorhabditis elegans
Lead Time	Made-to-order (14-16 weeks)
Gene Names	ZK1290.11