



atpF2 Antibody

Product Code	CSB-PA653242XA01FPY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q31RF4
Immunogen	Recombinant Synechococcus elongatus (strain PCC 7942) (Anacystis nidulans R2) atpF2 protein
Raised In	Rabbit
Species Reactivity	Synechococcus elongatus (strain PCC 7942) (Anacystis nidulans R2)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Synechococcus elongatus (strain PCC 7942) (Anacystis nidulans R2)
Lead Time	Made-to-order (14-16 weeks)
Target Names	atpF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.