



hemL Antibody

Product Code	CSB-PA657428XA01FPY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q31QJ2
Immunogen	Recombinant Synechococcus elongatus (strain PCC 7942) (Anacystis nidulans R2) hemL protein
Raised In	Rabbit
Species Reactivity	Synechococcus elongatus (strain PCC 7942) (Anacystis nidulans R2)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Synechococcus elongatus (strain PCC 7942) (Anacystis nidulans R2)
Lead Time	Made-to-order (14-16 weeks)
Target Names	hemL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.