



panC/cmK Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA653933XA01FPY |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q31P38 |
| Immunogen | Recombinant <i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2) panC/cmK protein |
| Raised In | Rabbit |
| Species Reactivity | <i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2) |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Product Type | Polyclonal Antibody |
| Immunogen Species | <i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2) |
| Lead Time | Made-to-order (14-16 weeks) |
| Target Names | panC/cmK |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |