



# murB Antibody

<b>Product Code</b>	CSB-PA310594XA01FPY
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P95837
<b>Immunogen</b>	Recombinant <i>Synechococcus elongatus</i> (strain PCC 7942) ( <i>Anacystis nidulans</i> R2) murB protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	<i>Synechococcus elongatus</i> (strain PCC 7942) ( <i>Anacystis nidulans</i> R2)
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	<i>Synechococcus elongatus</i> (strain PCC 7942) ( <i>Anacystis nidulans</i> R2)
<b>Lead Time</b>	Made-to-order (14-16 weeks)
<b>Target Names</b>	murB
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.