



ndhD1 Antibody

Product Code	CSB-PA653228XA01FPY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q31NA0
Immunogen	Recombinant Synechococcus elongatus (strain PCC 7942) (Anacystis nidulans R2) ndhD1 protein
Raised In	Rabbit
Species Reactivity	Synechococcus elongatus (strain PCC 7942) (Anacystis nidulans R2)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Synechococcus elongatus (strain PCC 7942) (Anacystis nidulans R2)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ndhD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.