



Recombinant *Xenopus laevis* Nuclear receptor subfamily 6 group A member 1 (nr6a1)

Product Code	CSB-YP016067XBE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P70033
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MDTWEDDRVD QRACLICGDR ATGLHYGIIS CEGCKGFFKR SICNKRVYRC SRDKNCVMSR KQRNRCQYCR LLKCLQMGMN RKAIREGMP GGRNKSIGPV QISDEEVERI MSGQEFEEEA NTSWSNNGDS DHSSPGNAVS ESNQPSPVST PSSSRSMELN GFSSLRDQYL GTPGPTHYQY LPHLFSFSAH PTIIPAQSRS LDPQSHTLIN QLLTAEDIEP LNTPLIEDG YKVTQSELFA LLCRLADELL FRQITWVKKL PFFCDLSIKD YTCLLSTTWQ ELILLSSITT YSKQIFGDLT DVTSKYSPE DELHRFSEDG MEVMERLIYL FRKFSQLKVS NEEYVCMKAI NFLNQDIQGI SSASQVEQLN KRYWYVCQDF TEYRYPHQPN RFPDLMMCLP EVRYIAGKLV NVPLEQLPLL FKAFLHSCKT IVTKE
Source	Yeast
Target Names	nr6a1
Protein Names	Recommended name: Nuclear receptor subfamily 6 group A member 1 Alternative name(s): Germ cell nuclear factor Short name= GCNF Short name= xGCNF
Expression Region	1-435
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.