



Amy58 Antibody

Product Code	CSB-PA515721XA01DKY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O18345
Immunogen	Recombinant Drosophila ananassae (Fruit fly) Amy58 protein
Raised In	Rabbit
Species Reactivity	Drosophila ananassae (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila ananassae (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Amy58
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.