



arm Antibody

Product Code	CSB-PA311606XA01MVC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02453
Immunogen	Recombinant Musca domestica (House fly) arm protein
Raised In	Rabbit
Species Reactivity	Musca domestica (House fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Musca domestica (House fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	arm
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.