



Amyrel Antibody

Product Code	CSB-PA529236XA01DMH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76459
Immunogen	Recombinant Drosophila serrata (Fruit fly) Amyrel protein
Raised In	Rabbit
Species Reactivity	Drosophila serrata (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila serrata (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Amyrel
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.