



ATPsynbeta Antibody

Product Code	CSB-PA632809XA01DMP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q24751
Immunogen	Recombinant Drosophila virilis (Fruit fly) ATPsynbeta protein
Raised In	Rabbit
Species Reactivity	Drosophila virilis (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila virilis (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ATPsynbeta
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.