



ORF4 Antibody

Product Code	CSB-PA631691XA01FQL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q27YG7
Immunogen	Recombinant Nora virus ORF4 protein
Raised In	Rabbit
Species Reactivity	Nora virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Nora virus
Lead Time	Made-to-order (14-16 weeks)
Target Names	ORF4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.