



CCHa2 Antibody

Product Code	CSB-PA122306XA01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8SXL2
Immunogen	Synthetic Peptide (27-39aa)
Species Reactivity	Drosophila melanogaster (Fruit fly)
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila melanogaster (Fruit fly)
Target Names	CCHa2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.