



GD21754 Antibody

Product Code	CSB-PA003313XA01DMJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B4Q9Z1
Immunogen	Recombinant Drosophila simulans (Fruit fly) GD21754 protein
Raised In	Rabbit
Species Reactivity	Drosophila simulans (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila simulans (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	GD21754
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.