



eIF3d1 Antibody

Product Code	CSB-PA463134XA01DKY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B3LY71
Immunogen	Recombinant Drosophila ananassae (Fruit fly) eIF3d1 protein
Raised In	Rabbit
Species Reactivity	Drosophila ananassae (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila ananassae (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	eIF3d1