



eIF3d2 Antibody

Product Code	CSB-PA455730XA01DMG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B4HHG8
Immunogen	Recombinant Drosophila sechellia (Fruit fly) eIF3d2 protein
Raised In	Rabbit
Species Reactivity	Drosophila sechellia (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila sechellia (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	eIF3d2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.