



EcR Antibody

Product Code	CSB-PA344227XA01CIP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49882
Immunogen	Recombinant Chironomus tentans (Midge) (Camptochironomus tentans) EcR protein
Raised In	Rabbit
Species Reactivity	Chironomus tentans (Midge) (Camptochironomus tentans)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Chironomus tentans (Midge) (Camptochironomus tentans)
Lead Time	Made-to-order (14-16 weeks)
Target Names	EcR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.