



ftsZ Antibody

Product Code	CSB-PA332608XA01WAR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45485
Immunogen	Recombinant Wolbachia sp. ftsZ protein
Raised In	Rabbit
Species Reactivity	Wolbachia sp.
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Wolbachia sp.
Lead Time	Made-to-order (14-16 weeks)
Target Names	ftsZ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.