




ftz-f1 Antibody

Product Code	CSB-PA339334XA01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33244
Immunogen	Recombinant Drosophila melanogaster (Fruit fly) ftz-f1 protein (790-1027aa)
Species Reactivity	Drosophila melanogaster (Fruit fly)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila melanogaster (Fruit fly)
Research Area	Others
Target Names	ftz-f1

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: ftz-f1 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 30 kDa

Observed band size: 30 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.