



FMRFa Antibody

Product Code	CSB-PA016032XA01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10552
Immunogen	Synthetic Peptide (204-208aa)
Species Reactivity	Drosophila melanogaster (Fruit fly)
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila melanogaster (Fruit fly)
Target Names	FMRFa
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.