



GM15071 Antibody

Product Code	CSB-PA453785XA01DMG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B4HXA6
Immunogen	Recombinant Drosophila sechellia (Fruit fly) GM15071 protein
Raised In	Rabbit
Species Reactivity	Drosophila sechellia (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila sechellia (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	GM15071
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.