



# M Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA965166XA01FQN   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | A7WNB1  |
| <b>Immunogen</b>           | Recombinant Drosophila melanogaster sigma virus (isolate Drosophila/USA/AP30/2005) (DMeISV) M protein |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Drosophila melanogaster sigma virus (isolate Drosophila/USA/AP30/2005) (DMeISV)                       |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)  |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | Non-conjugated  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                      |
| <b>Purification Method</b> | Antigen Affinity Purified   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Product Type</b>        | Polyclonal Antibody   |
| <b>Immunogen Species</b>   | Drosophila melanogaster sigma virus (isolate Drosophila/USA/AP30/2005) (DMeISV)                       |
| <b>Lead Time</b>           | Made-to-order (14-16 weeks)   |
| <b>Target Names</b>        | M   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                           |