



ix Antibody

Product Code	CSB-PA686323XA01BTT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5ISW3
Immunogen	Recombinant Bombyx mori (Silk moth) ix protein
Raised In	Rabbit
Species Reactivity	Bombyx mori (Silk moth)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	BOMMO(Silk moth)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ix
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.