



Lsp1gamma Antibody

Product Code	CSB-PA318191XA01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11997
Immunogen	Recombinant Drosophila melanogaster (Fruit fly) Lsp1gamma protein
Raised In	Rabbit
Species Reactivity	Drosophila melanogaster (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila melanogaster (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Lsp1gamma
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.