



Hsp83 Antibody

Product Code	CSB-PA519044XA01DLA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O02192
Immunogen	Recombinant Drosophila auraria (Fruit fly) Hsp83 protein
Raised In	Rabbit
Species Reactivity	Drosophila auraria (Fruit fly)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila auraria (Fruit fly)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Hsp83
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.